
GSK-3117391

产品编号: D51112

CAS: 1018673-42-1

分子式: C₂₂H₃₃N₃O₄

纯度: ≥98%

InChi: InChI=1S/C₂₂H₃₃N₃O₄/c26-20(25-28)13-12-18-11-10-16(14-23-18)15-24-21(17-6-2-1-3-7-17)2
2(27)29-19-8-4-5-9-19/h10-11,14,17,19,21,24,28H,1-9,12-13,15H2,(H,25,26)/t21-/m0/s1

InChi Key: AFDPFLDWOXXHQM-NRFXNRHFSA-N

Smiles: ONC(=O)CCC1=CC=C(CN[C@@H](C2CCCCC2)C(=O)OC2CCCCC2)C=N1

外观: 固体粉末

作用通路: HDAC

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: GSK-3117391, also known as GSK3117391A and HDAC-IN-3, is a potent HDAC inhibitor for treating chronic inflammatory disorders.