
SF-1670

产品编号: D51127

CAS: 345630-40-2

分子式: C₁₉H₁₇NO₃

纯度: ≥98%

InChi: InChI=1S/C₁₉H₁₇NO₃/c1-19(2,3)18(23)20-11-8-9-13-12-6-4-5-7-14(12)16(21)17(22)15(13)10-1
1/h4-10H,1-3H3,(H,20,23)

InChi Key: VZQDDSYKVYARDW-UHFFFAOYSA-N

Smiles: CC(C)(C)C(=O)NC1=CC2=C(C=C1)C1=CC=CC=C1C(=O)C2=O

外观: 固体粉末

作用通路: Autophagy

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: SF-1670 is a specific phosphatase and tensin homolog deleted on chromosome 10 (PTEN) inhibitor.