
TC-O 9311

产品编号: D51138

CAS: 444932-31-4

分子式: C₂₀H₁₉N₃O₄

纯度: ≥98%

InChi: InChI=1S/C₂₀H₁₉N₃O₄/c1-26-15-10-14(11-16(12-15)27-2)19(24)22-23-20(25)21-18-9-5-7-13-6
-3-4-8-17(13)18/h3-12H,1-2H3,(H,22,24)(H2,21,23,25)

InChi Key: KPTMSQHTGZMEFU-UHFFFAOYSA-N

Smiles: COC1C=C(C=C(C=1)OC)C(=O)NNC(=O)NC1=CC=CC2=CC=CC=C21

外观: 固体粉末

作用通路: GPR139

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: TC-O9311 is a potent GPR139 agonist (EC₅₀ = 39 nM in CHO-K1 cells expressing human GPR139).