
LTI-291

产品编号: D51170

CAS: 1919820-28-2

分子式: C₂₀H₃₀N₄O₂

纯度: ≥98%

InChi: InChI=1S/C₂₀H₃₀N₄O₂/c1-4-5-6-11-26-17-9-7-16(8-10-17)23-20(25)18-13-21-24-15(3)12-14(2)
22-19(18)24/h12-13,16-17H,4-11H₂,1-3H₃, (H,23,25)/t16-,17-

InChi Key: HZILSILAELSWKN-QAQDUYKDSA-N

Smiles: CCCCCO[C@@H]1CC[C@H](CC1)NC(=O)C1C=NN2C=1N=C(C)C=C2C

外观: 固体粉末

作用通路: Glucosidase

溶解性: DMSO : 50 mg/mL (139.48 mM; Need ultrasonic)

保存条件: Store in dry, dark place for one year.

产品介绍: LTI-291 is a glucocerebrosidase (Gcase) activator for treatment of Parkinson's disease (PD).