
DMAT

产品编号: D51193

CAS: 749234-11-5

分子式: C₉H₇Br₄N₃

纯度: ≥98%

InChi: InChi=1S/C₉H₇Br₄N₃/c1-16(2)9-14-7-5(12)3(10)4(11)6(13)8(7)15-9/h1-2H₃, (H,14,15)

InChi Key: SLPJGDQJLTYWCI-UHFFFAOYSA-N

Smiles: CN(C)C1NC2=C(N=1)C(Br)=C(Br)C(Br)=C2Br

外观: 固体粉末

作用通路: Casein Kinase

溶解性: DMSO : 50 mg/mL (104.87 mM; Need ultrasonic)

保存条件: Store in dry, dark place for one year.

产品介绍: DMAT is a potent and specific CK2 inhibitor with an IC₅₀ value of 130 nM.