
TM5275 sodium salt

产品编号: D51210

CAS: 1103926-82-4

分子式: C₂₈H₂₇ClN₃NaO₅

纯度: ≥98%

InChi: InChI=1S/C₂₈H₂₈ClN₃O₅.Na/c29-22-11-12-24(23(17-22)28(35)36)30-25(33)18-37-19-26(34)31-13-15-32(16-14-31)27(20-7-3-1-4-8-20)21-9-5-2-6-10-21;/h1-12,17,27H,13-16,18-19H₂, (H,30,33)(H,35,36);/q;+1/p-1

InChi Key: JSHSGBIWNPQCQZ-UHFFFAOYSA-M

Smiles: [Na+].[O-]
]C(=O)C1C=C(Cl)C=CC=1NC(=O)COCC(=O)N1CCN(CC1)C(C1C=CC=CC=1)C1C=CC=CC=1

外观: 固体粉末

作用通路: PAI-1

溶解性: Soluble in DMSO, not in water

保存条件: Store in dry, dark place for one year.

产品介绍: TM 5275 sodium salt is an orally bioavailable plasminogen activator inhibitor-1 (PAI-1) inhibitor.