

---

## WHI-P180

产品编号: D51235

CAS: 211555-08-7

分子式: C<sub>16</sub>H<sub>15</sub>N<sub>3</sub>O<sub>3</sub>

纯度: ≥98%

InChi: InChI=1S/C<sub>16</sub>H<sub>15</sub>N<sub>3</sub>O<sub>3</sub>/c1-21-14-7-12-13(8-15(14)22-2)17-9-18-16(12)19-10-4-3-5-11(20)6-10/h3-9,20H,1-2H3,(H,17,18,19)

InChi Key: BNDYIYYKEIXHNK-UHFFFAOYSA-N

Smiles: COC1C=C2C(NC3C=C(O)C=CC=3)=NC=NC2=CC=1OC

外观: 固体粉末

作用通路: EGFR

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: WHI-P180 is a potent EGFR and Cdk2 inhibitors with IC<sub>50</sub> of 4.0 and 1.0 μM, respectively.