

---

## Droxinostat

产品编号: D51261

CAS: 99873-43-5

分子式: C<sub>11</sub>H<sub>14</sub>ClNO<sub>3</sub>

纯度: ≥98%

InChi: InChi=1S/C11H14ClNO3/c1-8-7-9(12)4-5-10(8)16-6-2-3-11(14)13-15/h4-5,7,15H,2-3,6H2,1H3,(H,13,14)

InChi Key: JHSXDAWGLCZYSM-UHFFFAOYSA-N

Smiles: CC1C=C(Cl)C=CC=1OCCCC(=O)NO

外观: 固体粉末

作用通路: Apoptosis

溶解性: Soluble in DMSO, not in water

保存条件: Store in dry, dark place for one year.

产品介绍: Droxinostat, also known as NS 41080, is a selective inhibitor of HDAC3, HDAC6, and HDAC8. Droxinostat shows comparable inhibition of HDAC6 and HDAC8 (IC<sub>50</sub> = 2.47 and 1.46 μmol/L, respectively).