
Dp44mT

产品编号: D51308

CAS: 152095-12-0

分子式: C₁₄H₁₅N₅S

纯度: ≥98%

InChi: InChI=1S/C₁₄H₁₅N₅S/c1-19(2)14(20)18-17-13(11-7-3-5-9-15-11)12-8-4-6-10-16-12/h3-10H,1-2H₃, (H,18,20)

InChi Key: XOBIGRNRXCAMJQ-UHFFFAOYSA-N

Smiles: CN(C)C(=S)NN=C(C1C=CC=CN=1)C1C=CC=CN=1

外观: 固体粉末

作用通路: Apoptosis

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Dp44mT is an iron chelator, possessing DNA-damaging activity mediated by top2a inhibition.