
TA 02

产品编号: D51322

CAS: 1784751-19-4

分子式: C₂₀H₁₃F₂N₃

纯度: ≥98%

InChi: InChI=1S/C₂₀H₁₃F₂N₃/c21-15-7-5-13(6-8-15)18-19(14-9-11-23-12-10-14)25-20(24-18)16-3-1-2-4-17(16)22/h1-12H,(H,24,25)

InChi Key: QIFJOFNVIVQRNJ-UHFFFAOYSA-N

Smiles: FC1C=CC(=CC=1)C1N=C(NC=1C1C=CN=CC=1)C1C=CC=CC=1F

外观: 固体粉末

作用通路: Autophagy

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: TA 02 is a p38 MAPK inhibitor with IC₅₀ of 20 nM.