
TA-01

产品编号: D51323

CAS: 1784751-18-3

分子式: C₂₀H₁₂F₃N₃

纯度: ≥98%

InChi: InChI=1S/C₂₀H₁₂F₃N₃/c21-14-5-2-1-4-13(14)19-18(12-8-10-24-11-9-12)25-20(26-19)17-15(22)
6-3-7-16(17)23/h1-11H,(H,25,26)

InChi Key: QUGPSRHXNSBSL-UHFFFAOYSA-N

Smiles: FC1C=CC=CC=1C1N=C(NC=1C1C=CN=CC=1)C1C(F)=CC=CC=1F

外观: 固体粉末

作用通路: Autophagy

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: TA-01 is an inducer of cardiomyocyte differentiation from human embryonic stem cells following mesoderm induction.