
DL-AP4

产品编号: D51334

CAS: 20263-07-4

分子式: C₄H₁₀NO₅P

纯度: ≥98%

InChi: InChI=1S/C₄H₁₀NO₅P/c5-3(4(6)7)1-2-11(8,9)10/h3H,1-2,5H2,(H,6,7)(H2,8,9,10)

InChi Key: DDOQBQRIEWHWBT-UHFFFAOYSA-N

Smiles: NC(CCP(O)(O)=O)C(O)=O

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: DL-AP4 is an NMDA glutamate receptor antagonist.