
Bisindolylmaleimide I (GF109203X)

产品编号: D51355

CAS: 133052-90-1

分子式: C₂₅H₂₄N₄O₂

纯度: ≥98%

InChi: InChI=1S/C₂₅H₂₄N₄O₂/c1-28(2)12-7-13-29-15-19(17-9-4-6-11-21(17)29)23-22(24(30)27-25(23)31)18-14-26-20-10-5-3-8-16(18)20/h3-6,8-11,14-15,26H,7,12-13H₂,1-2H₃, (H,27,30,31)

InChi Key: QMGUOJYZJKLOLH-UHFFFAOYSA-N

Smiles: CN(C)CCCN1C=C(C2=C(C3=CNC4=CC=CC=C34)C(=O)NC2=O)C2=CC=CC=C12

外观: 固体粉末

作用通路: PKC

溶解性: Soluble in DMSO, not in water

保存条件: Store in dry, dark place for one year.

产品介绍: bisindolylmaleimide I is a potent and highly selective PCK inhibitor.