
p32-IN-M36

产品编号: D51380

CAS: 802555-85-7

分子式: C₂₃H₂₈N₈O₂

纯度: ≥98%

InChi: InChI=1S/C₂₃H₂₈N₈O₂/c32-22(30-18-6-1-4-16(14-18)20-24-10-11-25-20)28-8-3-9-29-23(33)31-19-7-2-5-17(15-19)21-26-12-13-27-21/h1-2,4-7,14-15H,3,8-13H₂, (H,24,25)(H,26,27)(H₂,28,30,32)(H₂,29,31,33)

InChi Key: HHFOXYBNULCICF-UHFFFAOYSA-N

Smiles: O=C(NC1=CC(=CC=C1)C1NCCN=1)NCCCNC(=O)NC1=CC(=CC=C1)C1NCCN=1

外观: 固体粉末

作用通路: PKC

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: p32-IN-M36 is an inhibitor of p32 mitochondrial protein, which binds directly to p32 and blocks p32 association with LyP-1.