
MDK-8650

产品编号: D51383

CAS: 736048-65-0

分子式: C₂₀H₂₅F₃N₆O₃

纯度: ≥98%

InChi: InChI=1S/C₂₀H₂₅F₃N₆O₃/c21-20(22,23)32-17-4-2-1-3-15(17)11-26-19-27-12-16(29(30)31)18(28-19)25-10-14-7-5-13(9-24)6-8-14/h1-4,12-14H,5-11,24H₂, (H₂,25,26,27,28)

InChi Key: HKOWATVSFKRXRW-UHFFFAOYSA-N

Smiles: NCC1CCC(CNC2=NC(NCC3C=CC=CC=3OC(F)(F)F)=NC=C2[N+](O-)=O)CC1

外观: 固体粉末

作用通路: PKC

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: MDK-8650, also known as PKC-θ inhibitor, is a selective PKC-θ inhibitor, with an IC₅₀ of 12 nM.