
KUN79359

产品编号: D51391

CAS: 257879-35-9

分子式: C₂₄H₂₁N₅O₂

纯度: ≥98%

InChi: InChI=1S/C₂₄H₂₁N₅O₂/c30-23-21(22(24(31)27-23)26-17-7-2-1-3-8-17)19-15-29(20-10-5-4-9-18(19)20)13-6-12-28-14-11-25-16-28/h1-5,7-11,14-16H,6,12-13H₂, (H₂,26,27,30,31)

InChi Key: KIWODJBCHRADND-UHFFFAOYSA-N

Smiles: O=C1NC(=O)C(NC2C=CC=CC=2)=C1C1=CN(CCCN2C=NC=C2)C2=CC=CC=C21

外观: 固体粉末

作用通路: Apoptosis

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: KUN79359, also known as RIP1 kinase inhibitor and PKCβ inhibitor, is a potent and selective inhibitor of PKCβ isozymes.