
NIBR3049

产品编号: D51402

CAS: 1260181-14-3

分子式: C₂₇H₂₅F₃N₄O₄

纯度: ≥98%

InChi: InChI=1S/C₂₇H₂₅F₃N₄O₄/c1-26(2,38)25(37)34-11-9-33(10-12-34)15-7-8-19(27(28,29)30)17(13-15)21-22(24(36)32-23(21)35)18-14-31-20-6-4-3-5-16(18)20/h3-8,13-14,31,38H,9-12H₂,1-2H₃,(H,32,35,36)

InChi Key: CLGRAWDGLMENOD-UHFFFAOYSA-N

Smiles: CC(C)(O)C(=O)N1CCN(CC1)C1=CC(C2=C(C3=CNC4=CC=CC=C43)C(=O)NC2=O)=C(C=C1)C(F)(F)F

外观: 固体粉末

作用通路: GSK-3

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: NIBR3049, also known as TCS-21311, is a potent and selective JAK3 inhibitor IC₅₀ values of 8 nM..