
Benclonidine

产品编号: D51535

CAS: 57647-79-7

分子式: C₁₆H₁₃Cl₂N₃O

纯度: ≥98%

InChi: InChi=1S/C₁₆H₁₃Cl₂N₃O/c17-12-7-4-8-13(18)14(12)20-16-19-9-10-21(16)15(22)11-5-2-1-3-6-11/h1-8H,9-10H₂, (H,19,20)

InChi Key: YHEMKMZUJKPOCO-UHFFFAOYSA-N

Smiles: O=C(C1C=CC=CC=1)N1CCN=C1NC1C(Cl)=CC=CC=1Cl

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Benclonidine is an Antihypertonic agent.