
HUN20688

产品编号: D51546

CAS: 246020-68-8

分子式: C₁₃H₁₀ClNO₅S

纯度: ≥98%

InChi: InChI=1S/C₁₃H₁₀ClNO₅S/c1-20-10-4-7-12(15(16)17)13(8-10)21(18,19)11-5-2-9(14)3-6-11/h2-8H,1H3

InChi Key: HURFJQIYDCUMPC-UHFFFAOYSA-N

Smiles: COC1=CC(=C(C=C1)[N+](=[O-])=O)S(=O)(=O)C1C=CC(Cl)=CC=1

外观: 固体粉末

作用通路: Others

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: HUN20688, also known as TXNRD1 inhibitor 1 ands TRi-1, is a potent, specific and irreversible inhibitor of cytosolic thioredoxin reductase 1 (TXNRD1).