
LYS-006

产品编号: D51569

CAS: 1799681-85-8

分子式: C₁₆H₁₄ClFN₆O₃

纯度: ≥98%

InChi: InChI=1S/C₁₆H₁₄ClFN₆O₃/c17-10-5-13(18)16(20-7-10)27-12-3-1-9(2-4-12)15-21-23-24(22-15)
8-11(19)6-14(25)26/h1-5,7,11H,6,8,19H₂, (H,25,26)/t11-/m0/s1

InChi Key: ZEGMEJVULDALSH-NSHDSACASA-N

Smiles: N[C@@H](CC(O)=O)CN1N=C(N=N1)C1C=CC(=CC=1)OC1=NC=C(Cl)C=C1F

外观: 固体粉末

作用通路: Aminopeptidase

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: LYS-006 is a Leukotriene A₄ hydrolase (LTA₄H) inhibitor for the treatment of Hidradenitis Suppurativa.