

---

## Benpyrine racemate

产品编号: D51585

CAS: 1333714-43-4

分子式: C<sub>16</sub>H<sub>16</sub>N<sub>6</sub>O

纯度: ≥98%

InChi: InChi=1S/C<sub>16</sub>H<sub>16</sub>N<sub>6</sub>O/c23-13-6-12(8-22(13)7-11-4-2-1-3-5-11)21-16-14-15(18-9-17-14)19-10-20-16/h1-5,9-10,12H,6-8H2,(H2,17,18,19,20,21)

InChi Key: HUWOMAVUXTXEKT-UHFFFAOYSA-N

Smiles: O=C1CC(CN1CC1C=CC=CC=1)NC1N=CN=C2NC=NC=12

外观: 固体粉末

作用通路: TNF Receptor

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: (Rac)-Benpyrine is a racemate of Benpyrine. (Rac)-Benpyrine is a potent and orally active TNF-α inhibitor and has the potential for TNF-α mediated inflammatory and autoimmune disease research.