
DSN28369

产品编号: D51605

CAS: 2413428-36-9

分子式: C₂₅H₃₁N₅O₈

纯度: ≥98%

InChi: InChI=1S/C₂₅H₃₁N₅O₈/c31-19(9-5-2-6-12-30-22(34)10-11-23(30)35)26-14-20(32)27-15-21(33)29-18(25(38)28-16-24(36)37)13-17-7-3-1-4-8-17/h1,3-4,7-8,10-11,18H,2,5-6,9,12-16H₂, (H,26,31)(H,27,32)(H,28,38)(H,29,33)(H,36,37)/t18-/m0/s1

InChi Key: DWPLKZYCOGLRPC-SFHVURJKSA-N

Smiles: OC(=O)CNC(=O)[C@H](CC1C=CC=CC=1)NC(=O)CNC(=O)CNC(=O)CCCCCN1C(=O)C=CC1=O

外观: 固体粉末

作用通路: ADC Linker

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: DSN28369 is a heterobifunctional linker, and useful to make antibody drug conjugate (ADC). DSN28369 contains a carboxy group, which can react with amino group or hydroxy group after activation. DSN28369 also contains a maleimido group, which can react with surfuryl group.