
DUN73423

产品编号: D51654

CAS: 2414373-42-3

分子式: C₁₉H₁₆N₆O

纯度: ≥98%

InChi: InChI=1S/C₁₉H₁₆N₆O/c20-18(26)15-3-1-14(2-4-15)16-12-23-25-10-7-17(24-19(16)25)22-11-13
-5-8-21-9-6-13/h1-10,12H,11H₂, (H₂,20,26)(H,22,24)

InChi Key: DICBAYNNPUCOOQ-UHFFFAOYSA-N

Smiles: NC(=O)C1C=CC(=CC=1)C1C=NN2C=CC(NCC3C=CN=CC=3)=NC2=1

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: DUN73423 is a RET/KDR-Selective Inhibitor of RETV804M Kinase with IC₅₀ = 7 nM. DUN73423 was listed in ACS Med Chem Lett . 2020 Feb 28;11(4):497-505 (compound 5).