
TLR4-C34

产品编号: D51677

CAS: 40592-88-9

分子式: C₁₇H₂₇NO₉

纯度: ≥98%

InChi: InChI=1S/C₁₇H₂₇NO₉/c1-8(2)24-17-14(18-9(3)19)16(26-12(6)22)15(25-11(5)21)13(27-17)7-23-10(4)20/h8,13-17H,7H2,1-6H3,(H,18,19)/t13-,14-,15-,16-,17+/m1/s1

InChi Key: KMIQMFHPUJUDMC-HHARLNAUSA-N

Smiles: CC(C)O[C@H]1O[C@H](COC(C)=O)[C@@H](OC(C)=O)[C@H](OC(C)=O)[C@H]1NC(C)=O

外观: 固体粉末

作用通路: Toll-like Receptor (TLR)

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: TLR4-C34 is a TLR4 inhibitor that inhibits TLR4 signaling in enterocytes and macrophages in vitro.