
Pomalidomide 4-PEG4-acid

产品编号: D51682

CAS: 2138440-81-8

分子式: C₂₄H₃₁N₃O₁₀

纯度: ≥98%

InChi: InChI=1S/C₂₄H₃₁N₃O₁₀/c28-19-5-4-18(22(31)26-19)27-23(32)16-2-1-3-17(21(16)24(27)33)25-7-9-35-11-13-37-15-14-36-12-10-34-8-6-20(29)30/h1-3,18,25H,4-15H₂, (H,29,30)(H,26,28,31)

InChi Key: CAANPUBZFFRWQP-UHFFFAOYSA-N

Smiles: OC(=O)CCOCCOCCOCCOCCNC1C=CC=C2C=1C(=O)N(C1CCC(=O)NC1=O)C2=O

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Pomalidomide 4-PEG4-acid is a functionalized cereblon ligand for PROTAC[®] research and development. It incorporates an E3 ligase ligand plus a PEG4 linker with terminal carboxylic acid ready for conjugation to a target protein ligand. It is part of a range of functionalized tool molecules for PROTAC R&D.