
NS 383

产品编号: D51688

CAS: 309711-59-9

分子式: C₁₉H₁₉N₃O₂

纯度: ≥98%

InChi: InChI=1S/C₁₉H₁₉N₃O₂/c1-2-22-9-8-13-14(12-6-4-3-5-7-12)10-15-17(16(13)11-22)20-19(23)18(15)21-24/h3-7,10,24H,2,8-9,11H₂,1H₃,(H,20,21,23)

InChi Key: CHVQAXYXYIDNDC-UHFFFAOYSA-N

Smiles: CCN1CC2C3NC(=O)/C(=N/O)/C=3C=C(C3C=CC=CC=3)C=2CC1

外观: 固体粉末

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: NS 383 is an ASIC blocker. Inhibition was observed at heteromeric ASIC channels. It attenuates pathophysiological nociceptive behaviors in CFA-inflamed and CCI rats.