
Thalidomide 4-ether-alkylC8-acid

产品编号: D51695

CAS: 2225148-51-4

分子式: C₂₂H₂₆N₂O₇

纯度: ≥98%

InChi: InChI=1S/C₂₂H₂₆N₂O₇/c25-17-12-11-15(20(28)23-17)24-21(29)14-8-7-9-16(19(14)22(24)30)3
1-13-6-4-2-1-3-5-10-18(26)27/h7-9,15H,1-6,10-13H₂, (H,26,27)(H,23,25,28)

InChi Key: JDVPWEYIVOQCDX-UHFFFAOYSA-N

Smiles: OC(=O)CCCCCCCCOC1C=CC=C2C=1C(=O)N(C1CCC(=O)NC1=O)C2=O

外观: 固体粉末

作用通路: Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Thalidomide 4-ether-alkylC8-acid is a cereblon ligand with alkyl linker and terminal acid for onward chemistry.