
Thalidomide 4-ether-alkylC4-acid

产品编号: D51696

CAS: 2169266-67-3

分子式: C₁₈H₁₈N₂O₇

纯度: ≥98%

InChi: InChI=1S/C₁₈H₁₈N₂O₇/c21-13-8-7-11(16(24)19-13)20-17(25)10-4-3-5-12(15(10)18(20)26)27-9-2-1-6-14(22)23/h3-5,11H,1-2,6-9H₂, (H,22,23)(H,19,21,24)

InChi Key: HCUNNTIALPNNPY-UHFFFAOYSA-N

Smiles: OC(=O)CCCCOC1C=CC=C2C=1C(=O)N(C1CCC(=O)NC1=O)C2=O

外观: 固体粉末

作用通路: E3 Ligase Ligand-Linker Conjugate

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: HUN 66673, also known as Thalidomide 4'-ether-alkylC4-acid is a cereblon ligand with alkyl linker and terminal acid for onward chemistry.