
Thalidomide 4-ether-alkylC3-acid

产品编号: D51700

CAS: 2169266-64-0

分子式: C₁₇H₁₆N₂O₇

纯度: ≥98%

InChi: InChI=1S/C₁₇H₁₆N₂O₇/c20-12-7-6-10(15(23)18-12)19-16(24)9-3-1-4-11(14(9)17(19)25)26-8-2-5-13(21)22/h1,3-4,10H,2,5-8H₂, (H,21,22)(H,18,20,23)

InChi Key: CQVOUORLTYFMHW-UHFFFAOYSA-N

Smiles: OC(=O)CCCOC1C=CC=C2C=1C(=O)N(C1CCC(=O)NC1=O)C2=O

外观: 固体粉末

作用通路: E3 Ligase Ligand-Linker Conjugate

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Thalidomide 4-ether-alkylC3-acid is a cereblon ligand with alkyl linker and terminal acid for onward chemistry.