

---

## Pomalidomide 4-PEG5-acid

产品编号: D51706

CAS: 2139348-63-1

分子式: C<sub>26</sub>H<sub>35</sub>N<sub>3</sub>O<sub>11</sub>

纯度: ≥98%

InChi: InChI=1S/C<sub>26</sub>H<sub>35</sub>N<sub>3</sub>O<sub>11</sub>/c30-21-5-4-20(24(33)28-21)29-25(34)18-2-1-3-19(23(18)26(29)35)27-7-9-37-11-13-39-15-17-40-16-14-38-12-10-36-8-6-22(31)32/h1-3,20,27H,4-17H<sub>2</sub>, (H,31,32)(H,28,30,33)

InChi Key: PMCQLKMOQUEEDR-UHFFFAOYSA-N

Smiles: OC(=O)CCOCCOCCOCCOCCOCCN<sub>1</sub>C=CC=C<sub>2</sub>C=1C(=O)N(C1CCC(=O)NC1=O)C<sub>2</sub>=O

外观: 固体粉末

作用通路: E3 Ligase Ligand-Linker Conjugate

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: Pomalidomide 4-PEG5-acid is a cereblon ligand with PEG linker and terminal acid for onward chemistry. This product has not formal name. For the convenience of scientific communication, we named it by combining its Inchi key (3 letters from the first letter of each section) with the last 5 digits of its CAS number or its molecule weight if its CAS number is not available.