
clAP1 ligand 1

产品编号: D51707

CAS: 2095244-42-9

分子式: C₃₁H₄₂N₄O₆S

纯度: ≥98%

InChi: InChI=1S/C₃₁H₄₂N₄O₆S/c1-19(34(5)30(40)41-31(2,3)4)27(38)33-25(20-11-7-6-8-12-20)29(39)
35-16-10-15-24(35)28-32-23(18-42-28)26(37)21-13-9-14-22(36)17-21/h9,13-14,17-20,24-25,36
H,6-8,10-12,15-16H₂,1-5H₃, (H,33,38)/t19-,24-,25-/m0/s1

InChi Key: SPXBMEKMIPOEMC-LQGLAIQGSA-N

Smiles: CN([C@@H](C)C(=O)N[C@@H](C1CCCCC1)C(=O)N1CCC[C@H]1C1=NC(=CS1)C(=O)C1=CC(O)=
CC=C1)C(=O)OC(C)(C)C

外观: 固体粉末

作用通路: Ligand for E3 Ligase

溶解性: Soluble in DMSO

保存条件: Store in dry, dark place for one year.

产品介绍: LCL 161, phenol is a functionalized IAP ligand for PROTACs.