

---

## Human Tau-4 Mouse mAb, HRP conjugated

Cat. No:	V5203-HRP
Name:	Human Tau-4 Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2A6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Buffer:	0.01M PBS (pH7.4).
Antigen:	Recombinant Human Tau-4 Protein