

bs-7550R

[Primary Antibody]

AUP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 550</p> <p>Target: AUP1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human AUP1/Ancient ubiquitous protein 1: 381-476/476.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Ancient ubiquitous protein 1 (AUP1) contains a domain with homology to the ancient conserved region of the archain 1 gene and a domain that may be involved in binding ubiquitin-conjugating enzymes. It has been shown to bind to the conserved membrane-proximal sequence of the cytoplasmic tail of integrin alpha(IIb) subunits.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9Y679</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 53 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	--	---