

---

## Recombinant E.coli BirA protein, N-His

Catalog Number: bs-42327P

Concentration: >1mg/ml

Species: Escherichia Coli

AA Seq: 1-321/321

Predicted MW: 37.4

Tags: N-His

Activity: Yes

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 50mM Imidazole (pH=6.8) with 50mM NaCl, 5% Glycerol, 5mM DTT.

Stored at -20°C.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.