
Recombinant Human Rhinovirus (HRV 3C) Protease , N-GST

Catalog Number: bs-42348P

Species: Human rhinovirus 14

AA Seq: 1-182/182

Predicted MW: 46.7 kDa

Tags: N-GST

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 150mM Tris-HCL (pH=8.0) with 150mM NaCL, 1mM DTT, 1mM EDTA, 50 % v/v glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.