

**bs-42330R**

**[ Primary Antibody ]**

## Mumps virus Non-structural protein V Rabbit pAb

**Bioss**  
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit


**Clonality:** Polyclonal

**Target:** Mumps virus Non-structural protein V

**Applications:** WB (1:500-2000)

**Reactivity:** Recombinant protein

### — VALIDATION IMAGES —



Sample: Lane 1: Recombinant Mumps virus Non-structural protein V(bs-42330P) Primary: Anti-Mumps virus Non-structural protein V (bs-42330R) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: kDa Observed band size: 36 kDa