
Recombinant mouse S100A7 protein

Catalog Number: bs-42023P

Concentration: >0.5mg/ml

Species: Mouse

AA Seq: 1-108/108

Predicted MW: 40.6 kDa

Tags: N-GST

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4) with 15mM GSH.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Fetal ear, skin, and tongue and human cell lines. Highly up-regulated in psoriatic epidermis. Also highly expressed in the urine of bladder squamous cell carcinoma (SCC) bearing patients.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.