
Recombinant mouse S100A7 protein

Catalog Number: bs-42021P

Concentration: >1mg/ml

Species: Mouse

AA Seq: 1-108/108

Predicted MW: 31 kDa

Tags: N-Trx-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCL (pH=8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Fetal ear, skin, and tongue and human cell lines. Highly up-regulated in psoriatic epidermis. Also highly expressed in the urine of bladder squamous cell carcinoma (SCC) bearing patients.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.