
Recombinant Immunoglobulin-degrading enzyme, N-His

Catalog Number: bs-42349P

Concentration: >1mg/ml

Species: Streptococcus pyogenes

AA Seq: 30-341/341

Predicted MW: 36.4

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.