

NF-L Mouse mAb

Cat. No: V5507
Name: NF-L Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 6B2
Host/Source: Mouse
Reactivity: (predicted: Human)
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human NF-L Protein ([bs-41217P](#))
Pairing Information:

ELISA		Detection									
		V5501	V5502	V5503	V5504	V5505	V5506	V5507	V5508	V5509	V5510
Capture	V5501	-	-	-	-	+	-	-	-	-	-
	V5502	-	-	-	-	-	-	-	-	-	-
	V5503	-	-	-	-	+	-	-	-	-	-
	V5504	-	+	-	-	-	-	-	-	-	-
	V5505	-	-	-	+	-	+	-	-	-	-
	V5506	-	-	-	-	-	-	-	-	-	-
	V5507	-	+	-	-	-	-	-	-	-	-
	V5508	-	-	-	-	-	-	-	-	+	-
	V5509	-	+++	-	-	-	-	-	-	-	-
	V5510	-	+++	-	-	-	-	-	-	-	-