
NF-L Mouse mAb

Cat. No:	V5508
Name:	NF-L Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3A5
Host/Source:	Mouse
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human NF-L Protein (bs-41217P)
Pairing Information:	
V5508- V5511& V5512	