

Alpha-synuclein Mouse mAb

Cat. No:	V5704
Name:	Alpha-synuclein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	14H8
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	Recombinant human Alpha-synuclein protein (bs-42005P)
Pairing Information:	

ELISA		Detection				
		V5702	V5703	V5704	V5705	V5706
Capture	V5702	-	-	-	-	-
	V5703	-	-	+	-	-
	V5704	+	-	-	-	-
	V5705	-	-	+++	-	+
	V5706	-	-	-	-	-