

Human UCHL1/PGP9.5 Mouse mAb

Cat. No: V5801
Name: Human UCHL1/PGP9.5 Mouse mAb
Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 11E4
Host/Source: Mouse
Subtype: IgG2a
Reactivity: Human
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: 0.01M PBS(pH7.4).
Applications: ELISA
Antigen: Recombinant Human PGP9.5 Protein ([bs-41257P](#))
Pairing Information:

ELISA		Detection				
		V5801	V5802	V5803	V5804	V5805
Capture	V5801	-	+	+	-	-
	V5802	+++	-	-	-	-
	V5803	+++	-	-	-	-
	V5804	-	+	-	-	-
	V5805	+	-	-	-	-