

**bs-39093R**

**[ Primary Antibody ]**

## Adenovirus Rabbit pAb

**Bioss**  
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Adenovirus

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml


**Storage:** 0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** WB (1:500-2000)

**Reactivity:** Adenovirus

### — VALIDATION IMAGES —



Sample: Lane 1: Adenovirus protein(bs-39093P)

Primary: Anti-Adenovirus (bs-39093R) at 1/1000

dilution Secondary: IRDye800CW Goat Anti-

Rabbit IgG at 1/20000 dilution Predicted band

size: kDa Observed band size: kDa