

**bsm-41591M**

**[ Primary Antibody ]**

## HBsAg Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG1

**Applications:** ELISA (1:5000-10000)

**Clonality:** Monoclonal

**Reactivity:** (predicted: HBsAg)

**Target:** HBsAg

**Immunogen:** Recombinant HBsAg protein: 2-119/400.

**Purification:** affinity purified by Protein A

**Predicted  
MW.:** 12/44 kDa

**Storage:** PBS, pH7.4.  
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated  
freeze/thaw cycles.

**Subcellular  
Location:** Cell membrane ,Cytoplasm

**Background:**