
PrV gB Rabbit pAb

Catalog Number: bs-20520R

Target Protein: PrV gB

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:PrV)

Predicted MW: 100 kDa

Entrez Gene: SuHV-1

Source: KLH conjugated synthetic peptide derived from human PrV gB: 820-913/913.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.