## bs-20520R

## [ Primary Antibody ]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## PrV gB Rabbit pAb

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: PrV gB

**Immunogen:** KLH conjugated synthetic peptide derived from PrV gB: 820-913/913.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Applications: ELISA** (1:5000-10000)

Reactivity: (predicted: PrV)

Predicted MW.: 100 kDa