
MMP9 Mouse mAb, Gold conjugated

Cat. No:	V5902-Gold
Name:	MMP9 Mouse mAb, Gold conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4B8
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Buffer:	0.01M PBS (pH7.4).
Antigen:	Recombinant human MMP9 protein