
Recombinant human Glyt1 protein

Catalog Number: bs-42369P

Concentration: >1mg/ml

Species: Human

AA Seq: 209-285/706

Predicted MW: 26.7 kDa

Tags: N-Trx-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 20mM Tris-HCl (pH=8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.