

## Recombinant mouse IL-2 protein, C-His

Catalog Number: bs-41224P

Concentration: >1mg/ml

Species: Mouse

AA Seq: 21-169/169

Predicted MW: 18 kDa

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC


Form: Liquid

Storage: 20mM Tris-HCl (pH=8.0) with 150mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** The protein encoded by this gene is a secreted cytokine that is important for the proliferation of T and B lymphocytes. The receptor of this cytokine is a heterotrimeric protein complex whose gamma chain is also shared by interleukin 4 (IL4) and interleukin 7 (IL7). The expression of this gene in mature thymocytes is monoallelic, which represents an unusual regulatory mode for controlling the precise expression of a single gene. The targeted disruption of a similar gene in mice leads to ulcerative colitis-like disease, which suggests an essential role of this gene in the immune response to antigenic stimuli.  
[provided by RefSeq, Jul 2008]

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.